

MOTOMED MUVI PASSIVE & ACTIVE TRAINER



SKU: FPM101100



PRODUCT DESCRIPTION

Purpose & Benefits

The **MOTomed Muvi Passive & Active Trainer** enables users to perform leg and arm training simultaneously, combining active and passive movement modes for comprehensive rehabilitation and exercise therapy. It allows independent adjustment of speed, resistance, and motor power for each limb group, promoting balanced and efficient therapy sessions for mobility, coordination, and muscle activation.

Features

- Simultaneous leg and arm training in both passive or active modes
- Independently adjustable movement speed and resistance for arms and legs
- Upgradeable via serial port to an FES Cycle using the HASOMED RehaMove2 Functional Electrical Stimulator
- Electronic foot insertion aid for ease of setup and use
- Passive speed range: 1 – 60 rpm; Active resistance: 0 – 20; Passive motor power: 1 – 10
- Therapy time preset adjustable from 0 to 120 minutes

- Movement protection, spasm control, and automatic safety stop for user safety
- Symmetry training promotes balanced limb engagement
- SmoothDriveSystem ensures quiet, gentle cycling with soft start and stop
- Therapy and motivation programs with detailed biofeedback during and after training
- Data export and import functions for personalised training programs
- Built-in games, videos, and slideshow functionality via USB for engagement
- Multi-language selection for user convenience

Specifications

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight
FPM101100	MOTOMED Muvi - Leg guides, Foot & Hand Fix	95 – 108 cm	60 cm	124 – 135 cm	135 kg	51 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au