



NEEKI WONDER SHEET PLUS - SHEET AND PILLOWCASE SET

\$242.00 - \$329.00 Price range: \$242.00 through \$329.00



SKU: WP-SETLS35G-H

PRODUCT DESCRIPTION

Purpose & Benefits

The **Wonder Sheet +PLUS by Neeki Sheet and Pillowcase Set** is a unique-fitted sheet design that increases mobility while reducing the effort and pain involved with turning over and getting in and out of bed. It is especially suited for those who require extra stability on the sides of the bed, such as residents in aged care facilities, where repositioning becomes easier and staff can reduce the risk of back injuries. The satin panel under the shoulders, back, and hips allows smoother sliding and pivoting, while the cotton/polyester borders provide grip and stability at the bed edges.

Features

- 'Glossy Charmeuse' satin panel for smooth movement
- Cotton/polyester borders for extra stability at the sides
- 50/50 cotton/polyester blend fabric
- Machine washable up to 71°C to meet aged care facility guidelines
- Pack includes 2 fitted sheets and 2 standard pillowcases
- Complies with AS/NZS4146:2000 Laundry Practice Guidelines

Specifications

- **Available Colours:** Green, Blue
- **Pack Sizes:** Includes 2 fitted sheets + 2 pillowcases
- **Long Single:** 91 x 203 x 35cm (Green, includes 2 pillowcases)
- **Queen:** 152 x 203 x 35cm (Blue or Green, includes 2 pillowcases)
- **King Single:** 107 x 203 x 35cm (Blue or Green, includes 2 pillowcases)
- **Mattress Fit:** 35cm depth fits mattresses up to 30cm; 45cm depth fits mattresses up to 40cm