



NEEN ANUFORM PROBE

\$112.00



SKU: PAT-091253640

PRODUCT DESCRIPTION

HOW IT WORKS

You can use the Anuform to stimulate your anal sphincter muscles and pelvic floor muscles, making them stronger, or to find out if you're contracting your muscles the way you should.

If it's to strengthen your muscles, use Anuform with Pericalm™. If it's feedback you're after, use your Anuform with the Peritone. The Anuform picks up electrical activity in the anal sphincter muscles and pelvic floor. This information then appears on your device, letting you know how you're doing with your exercises.

KEY FEATURES

- It's specially contoured for comfort
- Easy to clean
- Stays in place when you move around
- Enhanced EMG biofeedback results
- Fast stimulation of your muscles