



NICE STRETCH DORSAL

\$55.00 - \$57.00 Price range: \$55.00 through \$57.00



SKU: NS-52060-H

PRODUCT DESCRIPTION

Purpose & Benefits of Nice Stretch Dorsal

The Nice Stretch Dorsal night splint is designed to treat plantar fasciitis, Achilles tendonitis, and other related heel and tendon disorders. Unlike traditional posterior night splints, the Nice Stretch® Dorsal features a front-facing (dorsal) design that provides a more comfortable, less bulky fit and promotes effective stretching of the plantar fascia and Achilles tendon. This positioning helps alleviate the intense morning pain commonly associated with these conditions by maintaining a gentle, therapeutic stretch throughout the night. The dorsal design also allows for better airflow, making it cooler and more comfortable to wear during sleep.

Features

- **Dorsal Design:** Places the support on the top of the foot and along the anterior shin, avoiding bulkiness and heat retention associated with posterior models.
- **Adjustable Tension Straps:** Enables precise control over the amount of stretch applied, allowing for customization to individual therapeutic needs.
- **Lightweight Construction:** Made from lightweight materials, enhancing comfort and reducing the feeling of being restrained during sleep.
- **Minimalist Footprint:** The compact design reduces interference with bedding and makes it easier to sleep comfortably through the night.
- **Breathable Fabric:** Features breathable fabrics that help keep the foot cool, minimizing sweating and discomfort.
- **Non-Skid Pad:** Includes a non-skid pad on the sole to enhance safety when walking in the splint, such as during nighttime bathroom visits.

Specifications

- **Material:** High-quality, breathable fabric coupled with durable plastics for structural support.
- **Size:** Available in multiple sizes to accommodate different foot dimensions.
- **Care Instructions:** Easy to maintain; the fabric can be spot-cleaned as needed, and the straps can be wiped down with a damp cloth.