



ALRICK NINO CHILDREN'S BED



SKU: 7-NINO

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Nino Series offers children the opportunity to explore their surroundings while also ensuring their safety and well-being. The clear Perspex integrates seamlessly into any space, offering the individual an unobstructed perspective of their nearby environment while simultaneously safeguarding against potential falls or harm. Moreover, it is available in a variety of powder-coated hues and can be customized with taller door elements or doors on either side.

FEATURES

- Doors one side (750mm high perspex sides)
- 125mm braking and directional castors
- Battery backup
- Underbed lighting
- Backlit handset
- White laminate

SPECIFICATIONS

Parameter	Measurement/Detailed Info
Overall Dimensions	2200 x 1014 mm
Sleeping Surface	2000 x 870 mm
Safe Working Load (SWL)	200kg
Net Bed Weight	140kg
Low Position	390mm
High Position	830mm
Underbed Clearance	100mm (at lowest point)
Back Angle (°)	70

Parameter	Measurement/Detailed Info
Knee Angle (°)	44
Optional Features	Powdercoat colours, Higher doors (1100mm), Doors both sides (750mm only), Central locking castors

BROCHURE

[NINO-Series-Bed-By-Caremed-Alrick-2](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au