



## **NOVIS AREACARE COMFORT PRESSURE SURFACE MATTRESS - KING SINGLE**



**SKU:** BSPS6001-KS-H

### **PRODUCT DESCRIPTION**

#### **Purpose & Benefits**

The **Novis AreaCare Comfort Pressure Surface Mattress - King Single** range offers high-quality pressure redistribution and comfort for individuals in need of clinical-grade support. Designed for aged care, home care, or hospital settings, these mattresses are constructed from premium foam with medical-grade covers, offering excellent support for individuals at risk of pressure injuries.

Available in three options—**Comfort**, **Premium (Firm)**, and **Optimum (Soft)**—the AreaCare range allows caregivers to select the best surface according to the user's comfort needs, postural requirements, and pressure care level.

#### **Features**

- Designed to reduce pressure and support skin integrity
- High-grade foam construction with varying firmness levels
- Vapour-permeable, waterproof medical cover
- Fire-retardant and antimicrobial properties
- Suitable for hospital, aged care, and home use

- Available in standard, firm, and soft variations

## Specifications (All Variants)

- Size: King Single
- Length: 1980 mm
- Width: 1050 mm
- Height: 150 mm
- Maximum User Weight: 220 kg
- Manufacturer: Novis Healthcare
- Model Types: Comfort, Premium (Firm), Optimum (Soft)

For further information go to [medimart.com.au](https://www.medimart.com.au) or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)