

## NOVIS BARIATRIC ALTERNATING PRESSURE MATTRESS (APM) - EXPANDABLE WIDTH



**SKU:** APMPC-B01S

## PRODUCT DESCRIPTION

### Purpose & Benefits

The **Novis Bariatric Alternating Pressure Mattress (APM) - Expandable Width** delivers dynamic, therapeutic support tailored to bariatric users at high or very high risk of pressure injuries. With reinforced side bolsters and an alternating air-cell system, the mattress promotes circulation through cyclical inflation and deflation—helping to reduce skin breakdown and distribute pressure evenly, while offering exceptional lateral support and stability.

### Features

- Adaptive alternating pressure system designed for heavier individuals up to 454 kg
- Automatic Fowler Boost engages when the bed is adjusted to maintain sacral pressure relief
- Reinforced side bolsters enhance lateral stability and support
- Expandable width (36", 42", 48") and length (80", 85", 90") to fit larger or profiling beds

### Specifications

- **Dimensions (L × W × H):** 2032-2286 mm × 914-1219 mm × 250 mm
- **Safe Working Load:** 40-454 kg

- **Pressure Cycle:** Alternating transverse air cells (1 in 2 cycle)
- **Certifications:** IEC 60601 and CE-marked