



NOVIS BETTERLIVING MATTRESS OVERLAY SYSTEM

\$1,250.00



SKU: APMBL-L01

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The BetterLiving Mattress Overlay alternating pressure surface is a cost effective solution for the prevention and management of pressure related injuries.

Introducing the BetterLiving alternating pressure surfaces – a cost effective solution for the prevention and management of pressure related injuries.

Harnessing the latest technical features for improved clinical care and ease of use, BetterLiving offers both a mattress replacement and a mattress overlay option. Ideally suited to home care or aged care environments, these products can be used for people at moderate risk to high risk of pressure injury.

Safety

No more concerns over correct pressure level settings with our AutoCair feature – automatic weight sensing technology that delivers fail safe patient comfort levels.

Simplicity

Designed for basic ease of use and minimal maintenance. Easy to set up, simple to operate and straightforward to clean and maintain.

Quality

Conform to all relevant Australian standards. Made using the best quality materials, every product must pass our stringent Quality Control inspection and regular testing protocols.

Overlay

- Functional overlay for pressure prevention
- Automatic pressure setting
- Used over existing mattress for basic support
- Easy to operate and maintain
- 40-150 kg bariatric

FEATURES

- Feature rich overlay for pressure prevention
- 1:2 cell alternation overlay
- AutoCair feature – automatic weight sensing technology that delivers fail-safe patient comfort levels
- Designed to fit over an existing mattress or safety underlay
- Quick twist CPR for rapid deflation
- Carry bag included

RESOURCES

[View the Novis BetterLiving Mattress Overlay Brochure.](#)

SPECIFICATIONS

Length	Width	Height	SWL
2000mm	880mm	130mm	40kg-150kg