



# NOVIS BETTERLIVING MATTRESS OVERLAY SYSTEM

\$1,250.00 Original price was: \$1,250.00.<u>\$1,125.00</u>Current price is: \$1,125.00.



**SKU:** APMBL-L01

#### PRODUCT DESCRIPTION

#### **PURPOSE & BENEFITS**

The BetterLiving Mattress Overlay alternating pressure surface is a cost effective solution for the prevention and management of pressure related injuries.

Introducing the BetterLiving alternating pressure surfaces – a cost effective solution for the prevention and management of pressure related injuries.

Harnessing the latest technical features for improved clinical care and ease of use, BetterLiving offers both a mattress replacement and a mattress overlay option. Ideally suited to home care or aged care environments, these products can be used for people at moderate risk to high risk of pressure injury.

# Safety

No more concerns over correct pressure level settings with our AutoCair feature – automatic weight sensing technology that delivers fail safe patient comfort levels.

#### Simplicity

Designed for basic ease of use and minimal maintenance. Easy to set up, simple to operate and straightforward to clean and maintain.

### Quality

Conform to all relevant Australian standards. Made using the best quality materials, every product must pass our stringent Quality Control inspection and regular testing protocols.

# Overlay



- Functional overlay for pressure prevention
- Automatic pressure setting
- Used over existing mattress for basic support
- Easy to operate and maintain
- 40-150 kg bariatric

#### **FEATURES**

- Feature rich overlay for pressure prevention
- 1:2 cell alternation overlay
- AutoCair feature automatic weight sensing technology that delivers fail-safe patient comfort levels
- Designed to fit over an existing mattress or safety underlay
- · Quick twist CPR for rapid deflation
- Carry bag included

# **RESOURCES**

View the Novis BetterLiving Mattress Overlay Brochure.

### **SPECIFICATIONS**

Length	Width	Height	SWL	
2000mm	880mm	130mm	40kg-150kg	