



NOVIS BETTERLIVING MATTRESS REPLACEMENT SYSTEM

\$1,890.00



SKU: APMBL-R01

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

BetterLiving alternating pressure surfaces is a cost-effective solution for the prevention and management of pressure-related injuries. Ideally suited to home care or aged care environments, these products can be used for people at risk to high risk of pressure injury. Providing one of the best comfort and stress relieving bed surfaces on the market.

KEY FEATURES

Safety

No more concerns over correct pressure level settings with our AutoCair feature - automatic weight sensing technology that delivers fail-safe patient comfort levels.

Simplicity

Designed for basic ease of use and minimal maintenance. Easy to set up, simple to operate and straightforward to clean and maintain.

Quality

Conform to all relevant Australian standards. Made using the best quality materials, every product must pass our stringent Quality Control inspection and regular testing protocols.

EXTENDED FEATURES

- Essential features for pressure injury prevention and treatment
- Automatic pressure setting
- Replace existing mattress for superior therapy

- Cell in cell for enhanced user safety and comfort
- Easy to operate and maintain
- 40-180 kg

Harnessing the latest technical features for improved clinical care and ease of use, BetterLiving offers both a mattress replacement and a mattress overlay option.

SPECIFICATIONS

Patient Risk Category	High Risk
Support Surface APM	BetterLiving APM
Mattress Recommendations	Designed for high-end, home and age care needs. Offering a cost-effective solution while providing optimum patient comfort and safety. Easy to set up with automatic weight sensing and comfort
Clinical Features	2 cell cycle therapy Full Mattress Replacement Automatic weight sensing with manual comfort control
SWL (kg)	40-180kg
Dimensions (L X W X H) in mm	2000 x 880 x 200 mm