


NOVIS CAIRMAX DUO MATTRESS WITH FOAM AIR SUPPORT



SKU: FAMCM

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

CairMax by Novis – the new range of hybrid pressure surfaces, delivering optimum comfort in the treatment and prevention of pressure injuries.

The CairMax Duo is a versatile foam air support surface with the ability to alternate between static or dynamic support, in response to the needs of the patient. Mattress functions as a static pressure relieving surface without the control unit, however when clinically required, can ‘step up’ to provide dynamic support with the control unit connected.

KEY FEATURES

Mattress functions as a static pressure relieving surface without the control unit, however when clinically required, it can step up to provide dynamic support with the control unit connected.

- **Versatility** - Mattress and control unit can be purchased as a system or separately. The compact control unit can be easily connected for immediate step-up to dynamic therapy when required, with no disruption to patient comfort.
- **Safety** - No concerns over correct pressure level settings with AutoCair automatic weight sensing technology, delivering fail-safe patient comfort and support.
- **Simplicity** - This plug and play system is very simple to use and requires minimum training to operate, making it ideal for community or aged care use.
- **Comfort** - Combining durable foam core with alternating air cells embedded with castellated foam to deliver comfort and support for patients up to 230 kg. The castellated foam articulates when the mattress is profiled, while static head area foam offers stable support for enhanced

sleep.

SPECIFICATIONS

Product Name	Product Code	Length (mm)	Width (mm)	Height (mm)
CairMax Duo Foam Air Support Surface Mattress - Single	FAMCM-R01	2000	880	180
CairMax Duo Foam Air Support Surface Mattress - King Single	FAMCM-R01K	2000	1050	180
CairMax Duo Foam Air Support Surface Mattress Only - Single	FAMCM-M01	2000	880	180
CairMax Duo Foam Air Support Surface Mattress Only - King Single	FAMCM-M01K	2000	1050	180

Product Name	CairMax Duo Foam Air Support Surface Mattress
Patient Risk Category	At-Risk - High Risk
Wound Stage	1-3
Waterlow Score	15-20
Braden Score	<12
Support Surface APM	Cairmax Duo
Mattress Recommendations	Designed for high end, home and age care needs. Offers a versatile 'Step Up' solution in response to the needs of the patient. Static mattresses can easily be converted to dynamic therapy by adding the control unit.
Clinical Features	Hybrid 2 cell cycle therapy Mattress Replacement Automatic weight sensing with manual comfort control
SWL (kg)	230kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au