

NOVIS CAIRMAX DUO MATTRESS WITH FOAM AIR SUPPORT



\$2,376.00 - \$2,862.00Price range: \$2,376.00 through \$2,862.00



SKU: FAMCM

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

CairMax by Novis – the new range of hybrid pressure surfaces, delivering optimum comfort in the treatment and prevention of pressure injuries.

The CairMax Duo is a versatile foam air support surface with the ability to alternate between static or dynamic support, in response to the needs of the patient. Mattress functions as a static pressure relieving surface without the control unit, however when clinically required, can 'step up' to provide dynamic support with the control unit connected.

KEY FEATURES

Mattress functions as a static pressure relieving surface without the control unit, however when clinically required, it can step up to provide dynamic support with the control unit connected.

- Versatility Mattress and control unit can be purchased as a system or separately. The compact control unit can be easily connected for immediate step-up to dynamic therapy when required, with no disruption to patient comfort.
- Safety No concerns over correct pressure level settings with AutoCair automatic weight sensing technology, delivering fail-safe patient comfort and support.
- Simplicity This plug and play system is very simple to use and requires minimum training to operate, making it ideal for community or aged care use.
- Comfort Combining durable foam core with alternating air cells embedded with castellated foam to deliver comfort and support for patients up to 230 kg. The castellated foam articulates when the mattress is profiled, while static head area foam offers stable support for enhanced



sleep.

SPECIFICATIONS

Product Name	Product Code	Length (mm)	Width (mm)	Height (mm)
CairMax Duo Foam Air Support Surface Mattress - Single	FAMCM-R01	2000	880	180
CairMax Duo Foam Air Support Surface Mattress - King Single	FAMCM-R01K	2000	1050	180
CairMax Duo Foam Air Support Surface Mattress Only - Single	FAMCM-M01	2000	880	180
CairMax Duo Foam Air Support Surface Mattress Only - King Single	FAMCM-M01K	2000	1050	180

Product Name CairMax Duo Foam Air Support Surface Mattress

Patient Risk Category At-Risk - High Risk

Wound Stage 1-3 Waterlow Score 15-20 Braden Score <12

Support Surface APM Cairmax Duo

Designed for high end, home and age care needs. Offers a

versatile 'Step Up' solution in response to the needs of the patient.

Mattress Recommendations Static mattresses can easily be converted to dynamic therapy by

adding the control unit. Hybrid 2 cell cycle therapy Mattress Replacement

Clinical Features Automatic weight sensing

with manual comfort control

230kg SWL (kg)