



NOVIS PROCAIR BARIATRIC ALTERNATING PRESSURE MATTRESS SYSTEM

\$7,400.00 - \$8,400.00Price range: \$7,400.00 through \$8,400.00



SKU: APMPC-B04SW-H

PRODUCT DESCRIPTION

Purpose & Benefits

The **ProCair Bariatric Alternating Pressure Mattress System** is designed to provide advanced pressure care and superior comfort for bariatric patients at high risk of pressure injuries. This dynamic mattress system offers continuous pressure redistribution through alternating air cells, promoting healthy blood circulation and reducing skin shear. Suitable for use in hospitals, aged care, and home environments, it ensures safe and effective pressure injury prevention for users requiring high-capacity support.

Features

- Alternating pressure system designed to prevent and manage pressure injuries in bariatric patients
- · High-capacity design supports heavier users while maintaining stability and comfort
- Multi-stretch, vapour-permeable cover reduces friction and moisture for improved skin protection
- Digital control unit with adjustable pressure settings to suit individual user needs
- Auto-firm function for caregiver support during repositioning or nursing care



- CPR quick-release function for emergency deflation
- Low-pressure and power-failure alarms ensure patient safety and system reliability
- Quiet pump operation for patient comfort and restful sleep

Specifications

| Specification | Details |
|-------------------|--|
| Model | ProCair Bariatric Alternating Pressure Mattress System |
| Safe Working Load | Up to 450 kg |
| Mattress Size | 120 cm × 200 cm × 25 cm |
| Cell Type | Alternating air cells with static and auto-firm modes |
| Cover Material | Multi-stretch, vapour-permeable PU cover |
| Control Unit | Digital pump with adjustable cycle times and pressure levels |
| Power Supply | AC 230V, 50Hz |
| Alarm System | Low pressure and power failure alarms |
| CPR Function | Quick-release valve for rapid deflation |