

NOVIS PROCAIR BARIATRIC ALTERNATING PRESSURE MATTRESS SYSTEM



SKU: APMPC-B04SW-H

PRODUCT DESCRIPTION

Purpose & Benefits

The **ProCair Bariatric Alternating Pressure Mattress System** is designed to provide advanced pressure care and superior comfort for bariatric patients at high risk of pressure injuries. This dynamic mattress system offers continuous pressure redistribution through alternating air cells, promoting healthy blood circulation and reducing skin shear. Suitable for use in hospitals, aged care, and home environments, it ensures safe and effective pressure injury prevention for users requiring high-capacity support.

Features

- Alternating pressure system designed to prevent and manage pressure injuries in bariatric patients
- High-capacity design supports heavier users while maintaining stability and comfort
- Multi-stretch, vapour-permeable cover reduces friction and moisture for improved skin protection
- Digital control unit with adjustable pressure settings to suit individual user needs
- Auto-firm function for caregiver support during repositioning or nursing care

- CPR quick-release function for emergency deflation
- Low-pressure and power-failure alarms ensure patient safety and system reliability
- Quiet pump operation for patient comfort and restful sleep

Specifications

Specification	Details
Model	ProCair Bariatric Alternating Pressure Mattress System
Safe Working Load	Up to 450 kg
Mattress Size	120 cm × 200 cm × 25 cm
Cell Type	Alternating air cells with static and auto-firm modes
Cover Material	Multi-stretch, vapour-permeable PU cover
Control Unit	Digital pump with adjustable cycle times and pressure levels
Power Supply	AC 230V, 50Hz
Alarm System	Low pressure and power failure alarms
CPR Function	Quick-release valve for rapid deflation

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au