

NOVIS PROCAIR PLUS MATTRESS REPLACEMENT SYSTEM



SKU: APMPC-R0

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Novis ProCair Plus is being introduced as the latest range of alternating pressure solutions for the prevention and treatment of pressure injuries. Our feature-rich products offer a complete range of alternating pressure mattress solutions to suit any environment.

For high-risk patients, the Novis ProCair Plus serves as our premium solution, providing optimal patient therapy and comfort. It incorporates fully automated operational features, advanced safety features including a backup battery, and advanced error and service indicators.

The Novis ProCair Plus Mattress Replacement System, specifically designed for high-risk patients, ensures optimal patient therapy and comfort. It includes HeelCair, which offers increased protection and therapy for vulnerable heels. It has an independent heel zone consisting of narrower cells that can be disconnected as needed. To prevent "bottoming out," alternating air cells contain a lower "safety cell" that remains fully inflated, ensuring greater surface stability and patient safety.

With AutoCair, concerns about correct pressure level settings are alleviated. This feature utilizes continuous automatic weight-sensing technology to adjust the pressure settings automatically, maintaining the client within a safe pressure range for effective therapy. It ensures fail-safe patient comfort and support. The Fowler Boost mode automatically adjusts mattress pressure to provide increased support under the sacral area when the backrest is raised. Additionally, an inbuilt BackUp Battery unit offers up to 5 hours of emergency backup power in case of a power failure or accidental disconnection.

Equipped with the CairAlert feature, the ProCair Plus Mattress Replacement System includes an

optional warning alarm signal to indicate when a patient has left the mattress or bed. It also has an inbuilt service indicator, providing worry-free maintenance.

Our feature-rich products harness the best in technical innovation, offering a complete range of alternating-pressure mattress solutions for high-risk patients. They ensure optimal patient therapy and comfort.

KEY FEATURES

- **AutoCare:** Eliminate concerns about correct pressure level settings with our AutoCair feature. It utilizes continuous automatic weight sensing technology, delivering fail-safe patient comfort levels.
- **SafetyCare:** Enjoy the ultimate protection from "bottoming out." Alternating air cells include a lower "safety cell" that remains fully inflated, ensuring greater surface stability and patient safety.
- **HeelCare:** Provide increased protection and therapy for vulnerable heels with an independent heel zone consisting of narrower cells that can be disconnected as required.
- **FailSafe:** Have peace of mind with the inbuilt BackUp Battery unit, providing up to 5 hours of emergency backup power in the event of a power failure or accidental disconnection.
- **Automatic Fowler Boost:** Experience automatic mattress pressure adjustment to provide increased support under the sacral area when the backrest is raised.
- **CareAlert:** Stay informed when your patient leaves the bed with our inbuilt CairAlert feature. An optional warning alarm will indicate that the patient has gotten off the mattress.

Harnessing the best in technical innovation, our feature-rich products deliver a complete range of alternating pressure mattress solution for high-risk patients with optimal patient therapy and comfort.

- Sealed mattress base for optimal infection control
- Enlarged CPR valve for the acute care environment
- The real-time automatic pressure setting function adjusts to weight changes
- The backup battery provides up to five hours of continuous therapy without mains power
- Detailed error messages for specific servicing
- Automatic Fowler Boost when the bed is profiled
- Narrow heel cells can be disconnected individually for specialised treatment
- Bed exit alarm to notify the carer that the user has left the mattress
- Comes with a pump.

SPECIFICATIONS

Product Name	ProCair Plus Mattress Replacement System
Patient Risk Category	Very High Risk
Dimensions	Single (Standard): 2000mm x 880mm x 200mm King Single: 2000mm x 1050mm x 200mm

Mattress Recommendations	Designed for hospital and high-end age care needs. Alternating self-setting therapy in real-time Ensuring optimum pressure and shear reduction during bed articulation. Designed to prevent and treat pressure injuries Battery backup that ensures continuous therapy when mains power is not connected. Bed exit alarm for patients with cognitive dysfunction that pose a falls risk
Clinical Features	2 cell cycle therapy in real-time Full Mattress Replacement 5 Hours Battery Backup Automatic weight sensing with Automatic Fowler Boost in real-time! Heel Zone - Disconnectable
SWL (kg)	40-220kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au