

## NOVIS PROCAIR PRIME ADVANCED MATTRESS REPLACEMENT SYSTEM



**SKU:** APMPC-R03A-H

### PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The **ProCair Prime Advanced Mattress Replacement System** delivers high-end alternating pressure care for patients at very high risk of pressure injuries. Designed for hospitals and advanced aged care, it combines automatic weight-sensing technology, heel protection zones, and emergency backup power for uninterrupted support and therapy. With innovative AutoCair™ and HeelCair™ technologies, it provides real-time pressure adjustment, enhanced heel care, and reliable patient safety.

### FEATURES

- **HeelCair Technology:** Independent heel zone with narrower, disconnectable cells for advanced heel protection.
- **AutoCair Technology:** Automatically adjusts pressure based on patient weight for consistent comfort and safety.
- **Safety Cell Design:** Lower safety cells remain inflated to prevent “bottoming out” and improve stability.
- **Fowler Boost Mode:** Automatically increases support under the sacral area when the backrest is raised.
- **CairAlert Function:** Optional alarm alerts caregivers if a patient exits the bed.
- **Backup Battery:** Provides up to 5 hours of emergency operation in case of power loss.
- **Enlarged CPR Valve:** Enables quick deflation for emergency response situations.
- **Sealed Mattress Base:** Promotes hygiene and infection control.
- **Service Indicator:** Displays error messages for easy maintenance.
- **Top Cover:** Latex-free and easy to clean for improved patient safety.

## SPECIFICATIONS

- **Patient Risk Category:** Very High Risk
- **Recommended User Weight:** 40 – 220 kg
- **Warranty:** 3 years
- **Compliance:** ARTG Number 289458

### Single (Sealed Base)

- **Length:** 2000 mm
- **Width:** 880 mm
- **Height:** 200 mm

### King Single (Sealed Base)

- **Length:** 2000 mm
- **Width:** 1050 mm
- **Height:** 200 mm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)