

NOVIS PROCAIR 8 PRO MATTRESS REPLACEMENT SYSTEM



SKU: PTAM8211

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Introducing the Novis ProCair 8 Pro Mattress Replacement System: Economical Patient Therapy and Unmatched Comfort

The Novis ProCair 8 Pro Mattress Replacement System is a top choice for facilities prioritizing patient therapy and comfort while maintaining a strong infection control strategy. Designed specifically for individuals in the high to very high-risk pressure injury category, this system offers exceptional support and protection.

With alternating air cells, the Novis ProCair 8 Pro Mattress Replacement System ensures the ultimate protection from "bottoming out." Each cell contains a lower "safety cell" that remains fully inflated, providing enhanced surface stability and patient safety.

To address the specific needs of patients, the Manual Fowler Boost mode allows for adjustable mattress pressure, offering increased support under the sacral area when the backrest is raised. Worries about correct pressure level settings are alleviated with the innovative AutoCair technology. This automatic weight-sensing feature adjusts the pressure settings to maintain a safe pressure range for effective therapy, ensuring fail-safe patient comfort and support.

Maintenance is hassle-free with the inbuilt service indicator, providing peace of mind and efficient upkeep of the Novis ProCair 8 Pro Mattress Replacement System. It is available in two sizes: Standard and King Single. Simply select the desired size from the dropdown menu when placing an order.

KEY FUNCTIONS

Harnessing the best in technical innovation, our feature-rich products deliver a complete range of alternating pressure mattress solutions to suit any environment.

- ActivCare - Benefit from enhanced dynamic therapy and increased stability with the comfort of 1 in 3 alternation.
- AcuteCare - Maximum infection control
- SafetyCare - The ultimate protection from "bottoming out", alternating air cells contain a lower "safety cell" which remains fully inflated to provide greater surface stability and patient safety.
- Side Bolsters - Static side bolsters increase patient safety and sense of stability, by reducing the risk of cell collapse at side edges.
- HeelCare - Increased protection and therapy for vulnerable heels with an independent heel zone made up of narrower cells. Designed for individualised therapy to the at-risk heel and lower leg area, while preventing heels from sinking between the cells.
- ClimateCare - Delivering the maximum warmth and comfort from your alternating pressure system, our ClimateCare cells are constructed from pure thermoplastic polyurethane (TPU) material for enhanced heat retention qualities.

SPECIFICATIONS

- Safe Working Load - 40kg to 250kg
- System Alternation - 1 in 3
- Pressure Care Risk - High to Very High Risk
- Overall Length - 2000mm (single) / 2000mm (king single)
- Overall Width - 885mm (single) / 1050mm (king single)
- Overall Height - 200mm (single) / 200mm (king single)
- Warranty - 3 yrs

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au