



## OAPL SHORT AIR WALKER MOON BOOT



**SKU:** OAS-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The OAPL Short Air Walker is designed to immobilise the foot to decrease pain and swelling, promoting faster healing following injuries or surgical procedures.

The OAPL Short Air Walker features an integrated pump and valve to fill an air compartment within the lining of the boot. Air fills space around the foot and ankle to provide light compression as well as arch support, increasing stability of the injury. The air compression is easily adjusted to achieve a custom fit that is able to accommodate changes in swelling throughout the healing process. This Walker features a rocker sole for added comfort and ease of ambulation, as well as a non-slip sole to protect from slips, providing additional security in each step. The lining is a moisture wicking soft lining to provide ultimate comfort and breathability.

Designed to immobilise the foot or offload areas of the foot following metatarsal fractures, plantar fascia tears, bunioneictomies, forefoot soft tissue injuries.

### SIZING

| Size    | Womens (US) | Mens (US) | Footbed Length (cm) |
|---------|-------------|-----------|---------------------|
| X Small | 4 - 5.5     | Up to 4   | 22.5                |
| Small   | 6 - 8       | 4.5 - 6.5 | 26                  |
| Medium  | 8.5 - 10.5  | 7 - 9     | 29                  |
| Large   | 11 - 12+    | 9.5 - 12  | 30.5                |
| X Large | 12.5 - 14   | 32.5      |                     |

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)