




OAPL SHORT AIR WALKER MOON BOOT



SKU: OAS-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The OAPL Short Air Walker is designed to immobilise the foot to decrease pain and swelling, promoting faster healing following injuries or surgical procedures.

The OAPL Short Air Walker features an integrated pump and valve to fill an air compartment within the lining of the boot. Air fills space around the foot and ankle to provide light compression as well as arch support, increasing stability of the injury. The air compression is easily adjusted to achieve a custom fit that is able to accommodate changes in swelling throughout the healing process. This Walker features a rocker sole for added comfort and ease of ambulation, as well as a non-slip sole to protect from slips, providing additional security in each step. The lining is a moisture wicking soft lining to provide ultimate comfort and breathability.

Designed to immobilise the foot or offload areas of the foot following metatarsal fractures, plantar fascia tears, bunioneictomies, forefoot soft tissue injuries.

SIZING

| Size | Womens (US) | Mens (US) | Footbed Length (cm) |
|---------|-------------|-----------|---------------------|
| X Small | 4 - 5.5 | Up to 4 | 22.5 |
| Small | 6 - 8 | 4.5 - 6.5 | 26 |
| Medium | 8.5 - 10.5 | 7 - 9 | 29 |
| Large | 11 - 12+ | 9.5 - 12 | 30.5 |
| X Large | 12.5 - 14 | 32.5 | |

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au