



OAPL TALL AIR WALKER MOON BOOT

\$124.00



SKU: OAT-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Designed to immobilise the foot and ankle following injuries with the benefit of pneumatic support. Air pressure provides compression and comfort and is easily adjusted to provide a better fit and accommodate volume changes.

The high frame and rocker sole provides stability and reduced ankle movement for stable foot and ankle fractures, sprains, soft tissue injuries and postoperative applications.

SIZING

Size	Womens (US)	Mens (US)	Footbed Length (cm)
X Small	4 - 5.5	Up to 4	22.5
Small	6 - 8	4.5 - 6.5	26
Medium	8.5 - 10.5	7 - 9	29
Large	11 - 12+	9.5 - 12	30.5
X Large	12.5 - 14		