



## OAPL TALL ANKLE WALKER MOON BOOT

**\$102.00**



**SKU: ORT-**

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The OAPL Standard Walker is designed to immobilise the foot and ankle to decrease pain and swelling, promoting faster healing following injuries or surgical procedures.

### FEATURES

- Patented metal interlocker uprights for increased support and durability of the boot. Designed to stabilise the ankle at 90° and reduce movement of the ankle.
- Velcro closures to allow custom fit to accommodate changes in volume after injury.
- A rocker sole for added comfort and ease of ambulation with a high grip sole to protect from slips.
- Breathable, moisture wicking, hand washable liner.
- EVA foot bed for impact absorption and comfort.

### SIZE

Size	Womens (US)	Mens (US)	Footbed Length (cm)
------	-------------	-----------	---------------------

X Small	4 - 5.5	Up to 4	22.5
Small	6 - 8	4.5 - 6.5	26
Medium	8.5 - 10.5	7 - 9	29
Large	11 - 12	9.5 - 12	30.5
X Large	12.5 +	12.5 - 14	32.5