




OAPL TENNIS ELBOW STRAP - MEDIUM (BLACK)

~~\$20.00~~ Original price was: \$20.00. **\$13.00** Current price is: \$13.00. Sale Available Online Only! While Stocks Last



SKU: 46008ME

PRODUCT DESCRIPTION

Purpose & Benefits

The **OAPL Tennis Elbow Strap** is designed to ease discomfort from repetitive strain injuries such as tennis or golfer's elbow. By restricting the full contraction of forearm muscles when making a fist, it helps reduce tendon load, alleviate pain, and support the healing process—all while allowing for continued activity.

Features

- Compression strap applying pressure to the tendon area to reduce strain
- Nylon support with slight elasticity for flexible yet secure wear
- Adjustable tension via Velcro strap for personalised fit

Specifications

- **Brand:** OAPL
- **Colour:** Black

- **Size:** Medium
- **Material:** Nylon with slight elastic properties
- **Fastener:** Adjustable Velcro
- **Care Instructions:** Hand-wash gently with mild soap; air-dry thoroughly; do not bleach, tumble-dry, or iron