




OAPL TENNIS ELBOW STRAP - MEDIUM (BLACK)



SKU: 46008ME

PRODUCT DESCRIPTION

Purpose & Benefits

The **OAPL Tennis Elbow Strap** is designed to ease discomfort from repetitive strain injuries such as tennis or golfer's elbow. By restricting the full contraction of forearm muscles when making a fist, it helps reduce tendon load, alleviate pain, and support the healing process—all while allowing for continued activity.

Features

- Compression strap applying pressure to the tendon area to reduce strain
- Nylon support with slight elasticity for flexible yet secure wear
- Adjustable tension via Velcro strap for personalised fit

Specifications

- **Brand:** OAPL
- **Colour:** Black

- **Size:** Medium
- **Material:** Nylon with slight elastic properties
- **Fastener:** Adjustable Velcro
- **Care Instructions:** Hand-wash gently with mild soap; air-dry thoroughly; do not bleach, tumble-dry, or iron

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au