



OPPO MEDICAL ARM SLING



SKU: OPP308730-H

PRODUCT DESCRIPTION

Purpose & Benefits

The **OPPO Arm Sling (Model 3087)** is designed to provide essential support and stabilisation for injured or immobilised arms. It features an adjustable foam shoulder pad that enhances comfort by distributing pressure evenly, and its ultralight, breathable material reduces discomfort and improves ventilation during recovery.

Features

- Adjustable foam pad for customised comfort and pressure relief
- Hook-and-loop fasteners for quick, secure and precise fitting
- Ultra-lightweight, breathable fabric for enhanced comfort
- Hand-washable in cold water with mild soap; air-dry away from heat and sunlight

Specifications

- **Brand:** OPPO Medical
- **Material Composition:** 80 % polyester, 10 % nylon, 10 % foam

- **Colour:** Blue
- **Care Instructions:** Hand wash in cold water with mild soap; air-dry away from heat and sunlight; do not bleach, iron, dry-clean or tumble dry

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au