



OPPO INDUSTRIAL BACK SUPPORT

~~\$55.00~~ Original price was: ~~\$55.00~~. ~~\$40.00~~ Current price is: ~~\$40.00~~. Sale Available Online Only! While Stocks Last



SKU: OPP216940-H

PRODUCT DESCRIPTION

Purpose & Benefits

The **OPPO Industrial Back Support** is expertly engineered to deliver robust spinal and lumbar support, especially suited for physically demanding settings such as manufacturing, logistics, or construction. Its strategically reinforced straps and lightweight ventilated elastic provide firm compression while maintaining breathability and flexibility. Four flexible stays offer essential structural support, helping to maintain posture and reduce strain on the lower back during extended periods of lifting, bending, or standing. With a 9-inch back pad that distributes pressure evenly across the lumbar region, this support reduces muscle fatigue, promotes correct spinal alignment, and boosts endurance and comfort throughout demanding workdays.

Features

- 9-inch back panel delivers extended coverage and balanced pressure across the lower back
- Four flexible stays enhance spinal support while allowing natural movement
- Reinforced straps ensure a secure and stable fit during active use
- Lightweight, ventilated elastic fabric offers firm compression without compromising breathability
- Adjustable shoulder straps simplify application and allow for customised fitting

Specifications

- **Brand:** OPPO Medical
- **Available Sizes:** S, M, L, XXL
- **Colour:** Black
- **Material Composition:** 70 % Nylon, 27 % Rubber, 3 % Plastic board (stays)
- **Care Instructions:** Hand-wash in cold water with mild soap; air-dry away from heat and sunlight; do not bleach, iron, dry-clean or tumble-dry