



## OPPO KNEE IMMOBILISER LENGTH 20IN



**SKU:** OPP4030-2050-H

## PRODUCT DESCRIPTION

### Purpose & Benefits

The **OPPO Knee Immobiliser** is expertly designed to deliver firm yet gentle stabilisation to the knee joint, particularly useful following injury, surgery, or in cases of instability and ligament strain. Its open-front design and adjustable upper and lower straps make it simple to fit and adapt, ensuring a comfortable and secure application. Embedded aluminium stays along the sides maintain alignment and minimise unwanted movement, while the patella relief cutout helps reduce pressure on the kneecap. Overall, this support provides confident, pain-reducing immobilisation that promotes healing and reduces strain during recovery.

### Features

- Five aluminium stays deliver structured support and effective immobilisation
- Upper and lower stabilising straps allow fine-tuned tensioning for optimal fit and comfort
- Open-front design enables easy fitting and adjustment, even with limited mobility
- Patella pressure-relief zone helps reduce discomfort around the kneecap

### Specifications

- **Brand:** OPPO Medical

- **Available Sizes:** S, M, L
- **Colour:** Beige
- **Material Composition:** 80 % polyester, 15 % foam, 5 % nylon
- **Care Instructions:** Remove stays before hand washing in cold water with mild soap; air-dry away from heat and sunlight; do not bleach, iron, dry-clean or tumble-dry

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)