



## **OPPO KNEE IMMOBILISER LENGTH 20IN**

\$57.00 Original price was: \$57.00.\$44.00 Current price is: \$44.00.Sale Available Online Only! While Stocks Last



SKU: OPP4030-2050-H

### PRODUCT DESCRIPTION

## **Purpose & Benefits**

The **OPPO Knee Immobiliser** is expertly designed to deliver firm yet gentle stabilisation to the knee joint, particularly useful following injury, surgery, or in cases of instability and ligament strain. Its open-front design and adjustable upper and lower straps make it simple to fit and adapt, ensuring a comfortable and secure application. Embedded aluminium stays along the sides maintain alignment and minimise unwanted movement, while the patella relief cutout helps reduce pressure on the kneecap. Overall, this support provides confident, pain-reducing immobilisation that promotes healing and reduces strain during recovery.

#### **Features**

- Five aluminium stays deliver structured support and effective immobilisation
- Upper and lower stabilising straps allow fine-tuned tensioning for optimal fit and comfort
- Open-front design enables easy fitting and adjustment, even with limited mobility
- Patella pressure-relief zone helps reduce discomfort around the kneecap

# **Specifications**

• Brand: OPPO Medical



• Available Sizes: S, M, L

• Colour: Beige

• Material Composition: 80 % polyester, 15 % foam, 5 % nylon

• Care Instructions: Remove stays before hand washing in cold water with mild soap; air-dry away from heat and sunlight; do not bleach, iron, dry-clean or tumble-dry