



OPPO TOTAL STABILITY ANKLE BRACE

~~\$64.00~~ Original price was: ~~\$64.00~~. ~~\$46.00~~ Current price is: ~~\$46.00~~. Sale Available Online Only! While Stocks Last



SKU: OPP400560-H

PRODUCT DESCRIPTION

Purpose & Benefits

The **OPPO Total Stability Ankle Brace** delivers exceptional support and protection for weak or injured ankles. Its innovative heel-lock strapping system combined with Figure-8 straps effectively restricts inversion and eversion movements that often lead to sprains. Removable plastic splints enhance customisable support during different stages of activity or recovery, while the universally adaptable design accommodates both left and right ankles with ease. Crafted for comfort, the soft padding at the heel and tongue, alongside its lace-up application, minimises discomfort and eliminates the need for complex taping—making it perfect for active lifestyles and reliable rehabilitation support.

Features

- Heel-lock strapping stabilises the ankle against inversion and eversion sprains
- Non-stretch Figure-8 straps enhance stability during movement
- Removable bilateral plastic splints allow adjustable support as needed
- Universal fit accommodates either ankle
- Lace-up design makes application quick and intuitive

- Padded heel and tongue deliver enhanced comfort during prolonged wear

Specifications

- **Brand:** OPPO Medical
- **Sizes Available:** X-, XL
- **Colour:** Typically beige or neutral tones (depending on stock)
- **Care Recommendations:** Recommended to be worn with a sock to improve comfort and minimise irritation