



OPPO TOTAL STABILITY ANKLE BRACE



SKU: OPP400560-H

PRODUCT DESCRIPTION

Purpose & Benefits

The **OPPO Total Stability Ankle Brace** delivers exceptional support and protection for weak or injured ankles. Its innovative heel-lock strapping system combined with Figure-8 straps effectively restricts inversion and eversion movements that often lead to sprains. Removable plastic splints enhance customisable support during different stages of activity or recovery, while the universally adaptable design accommodates both left and right ankles with ease. Crafted for comfort, the soft padding at the heel and tongue, alongside its lace-up application, minimises discomfort and eliminates the need for complex taping—making it perfect for active lifestyles and reliable rehabilitation support.

Features

- Heel-lock strapping stabilises the ankle against inversion and eversion sprains
- Non-stretch Figure-8 straps enhance stability during movement
- Removable bilateral plastic splints allow adjustable support as needed
- Universal fit accommodates either ankle
- Lace-up design makes application quick and intuitive

- Padded heel and tongue deliver enhanced comfort during prolonged wear

Specifications

- **Brand:** OPPO Medical
- **Sizes Available:** X-, XL
- **Colour:** Typically beige or neutral tones (depending on stock)
- **Care Recommendations:** Recommended to be worn with a sock to improve comfort and minimise irritation

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au

