




OSCAR BARWON RECLINER LIFT CHAIR - SINGLE MOTOR

\$2,729.00 - \$2,857.00



SKU: 31056-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The new Oscar Barwon is part of the Australian-made Oscar deluxe lift chair range.

The chair features a velcro-attached three-roll back rest, offering excellent support to lower, middle and upper back. The top two rolls are filled with a fibrous filling that has a wonderful soft feel, whilst the bottom roll is a high quality foam that provides excellent lower lumbar support. The backrest is fully adjustable to meet individual needs of comfort, and cater to conditions including Scoliosis, Kyphosis and Arthritis.

Seat dimensions can be customised to suit individual requirements in addition to a range of optional extras including massage and pressure care seating.

Available in standard Single or Dual Motor in Single or Dual Motor configuration in a range of different fabrics to suit any home.

Available with or without show wood on arms to assist rising and protect the fabric or as a waterfall arm.

This model is the single motor version.

SIZES

- **Mini:** Suitable for users from approximately 4'10"- 5' 2" (147-157cm) with a weight capacity up to 130kg.
- **A - Small:** Suitable for users from approximately 5' 2"- 5' 6" (157-168cm) with a weight

capacity up to 130kg.

- **B - Medium:** Suitable for users from approximately 5'4"-5'10" (163-177cm) with a weight capacity up to 130kg.
- **C - Large:** Suitable for users from approximately 5'11"-6'5" (180-195cm) with a weight capacity up to 130kg.

OPTIONS & MODIFICATIONS

- Seat depth
- Seat height
- Seat width
- Armrest style
- Castor size (can also choose to remove castors)
- Armrest height
- Backrest height
- Footrest length
- Battery Backup
- Heat & Massage options
- Cushioning options (including ROHO cushion cutout)
- Rear push handle
- Gluteal recess
- Disable lift function
- Controller side

SPECIFICATIONS

- Barwon Zero M Lift recliner - Mini :
 - Seat height: 456mm
 - Seat Depth: 445mm
 - Seat Width: 540mm
 - Back Height: 636mm
 - Scap Height: 429mm
 - User Weight Capacity: 130kg
- Barwon A Lift recliner:
 - Seat height: 476mm
 - Seat Depth: 490mm
 - Seat Width: 540mm
 - Back Height: 647mm
 - Scap Height: 440mm
 - User Weight Capacity: 130kg
- Barwon B Lift recliner:
 - Seat height: 496mm
 - Seat Depth: 520mm

- Seat Width: 540mm
- Back Height: 688mm
- Scap Height: 481mm
- User Weight Capacity: 130kg
- Barwon C Lift recliner:
 - Seat height: 546mm
 - Seat Depth: 545mm
 - Seat Width: 540mm
 - Back Height: 718mm
 - Scap Height: 511mm
 - User Weight Capacity: 130kg

WARRANTY

2 Years Warranty: Mechanism, fabric and electronics

5 Years Warranty: Timber Frame and Foam.