



## OSCAR BARWON RECLINER LIFT CHAIR

**\$2,470.00 - \$3,986.00** Price range: \$2,470.00 through \$3,986.00



**SKU:** 31056M-2M-MA-CA-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Oscar Barwon Lift Recliner is a comfortable, Australian-made chair designed to provide ergonomic support and customisable sizing for users with different postural needs. It features a soft, three-roll backrest for targeted support of the lower, middle and upper back, with the lower roll constructed of firm foam for lumbar reinforcement.

Customisation options are available to suit a wide range of needs including seat dimensions, cushioning, and arm styles. Available in multiple fabric or vinyl finishes and with or without show wood armrests for aesthetic and practical enhancement.

This model is available in the following configurations:

- Single Motor
- Single Motor with Comfort Pack
- Dual Motor
- Dual Motor with Comfort Pack

### SIZES

- **Mini:** Suitable for users from approximately 4'10"-5'2" (147-157cm), up to 130kg.
- **Size A:** Suitable for users from approximately 5'2"-5'6" (157-168cm), up to 130kg.
- **Size B:** Suitable for users from approximately 5'4"-5'10" (163-177cm), up to 130kg.

- **Size C:** Suitable for users from approximately 5'11"-6'5" (180-195cm), up to 130kg.

## OPTIONS & MODIFICATIONS

- Seat depth
- Seat height
- Seat width
- Armrest style
- Castor size (or remove castors)
- Armrest height
- Backrest height
- Footrest length
- Battery backup
- Heat & massage
- Cushioning (incl. ROHO cutout)
- Rear push handle
- Gluteal recess
- Disable lift function
- Controller side

## SPECIFICATIONS

- Barwon Mini Lift Recliner:
  - Length: 900mm
  - Width: 830mm
  - Height: 1120mm
  - User Weight Capacity: 130kg
- Barwon A Lift Recliner:
  - Length: 900mm
  - Width: 830mm
  - Height: 1120mm
  - User Weight Capacity: 130kg
- Barwon B Lift Recliner:
  - Length: 930mm
  - Width: 830mm
  - Height: 1160mm
  - User Weight Capacity: 130kg
- Barwon C Lift Recliner:
  - Length: 950mm
  - Width: 830mm
  - Height: 1220mm
  - User Weight Capacity: 130kg

## **WARRANTY**

2 Years Warranty: Mechanism, fabric and electronics

5 Years Warranty: Timber Frame and Foam.