



OSCAR EASY GLIDE DINING CHAIR

\$762.00 Original price was: \$762.00.\$686.00 Current price is: \$686.00.



SKU: MAEA-DT

PRODUCT DESCRIPTION

OSCAR EASY GLIDE

Dining is made easy for everyone. Lightweight with a simple weight-activated gliding action, this Oscar Easy Glide dining chair increases independence and reduces injuries. Whether you're in a commercial or domestic setting you can't go past the Easy Glide.

FEATURES

- Allows the user or carer to effortlessly roll into the table.
- 50mm Rear self-locking wheels (release when seated locks when you stand for stability).
- Light and Rigid Aluminium Frame.
- 38mm front fixed wheels.
- Easy clean wood grain finish.
- 120kg weight capacity.
- High <u>resilience foam</u>.
- Available with matching arm pads.

SPECIFICATIONS

Seat Height: 480mmSeat Width: 500mmSeat Depth: 450mmChair Height: 830mm



Back Height: 370mmChair Width: 570mmArm Height: 690mm

WARRANTY

2 Years Warranty: Mechanism, fabric and electronics

5 Years Warranty: Timber Frame and Foam.

Have a look at other <u>utility chairs</u>.