



OSCAR EASY GLIDE DINING CHAIR

\$762.00



SKU: MAEA-DT

PRODUCT DESCRIPTION

OSCAR EASY GLIDE

Dining is made easy for everyone. Lightweight with a simple weight-activated gliding action, this Oscar Easy Glide dining chair increases independence and reduces injuries. Whether you're in a commercial or domestic setting you can't go past the Easy Glide.

FEATURES

- Allows the user or carer to effortlessly roll into the table.
- 50mm Rear self-locking wheels (release when seated locks when you stand for stability).
- Light and Rigid Aluminium Frame.
- 38mm front fixed wheels.
- Easy clean wood grain finish.
- 120kg weight capacity.
- High [resilience foam](#).
- Available with matching arm pads.

SPECIFICATIONS

- Seat Height: 480mm
- Seat Width: 500mm
- Seat Depth: 450mm
- Chair Height: 830mm

- Back Height: 370mm
- Chair Width: 570mm
- Arm Height: 690mm

WARRANTY

2 Years Warranty: Mechanism, fabric and electronics

5 Years Warranty: Timber Frame and Foam.

Have a look at other [utility chairs](#).