



OSCAR SLIMLINE RECLINER LIFT CHAIR - 4 MOTORS



SKU: SLM-HL-04-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Slimline lift recliner has a modern appearance. As well as having all the features of a lift recliner, it has an adjustable headrest and adjustable lumbar support.

Another function of the Slimline is a zero-gravity mechanism. A zero-gravity chair allows the user to get into a comfortable position, relaxing the spine, improving blood flow and lung capacity.

The lift assist function tilts and lifts the chair forward to aid with standing and relieve the strain of standing from a low position.

Showwood arms offer further support, allowing you to hold onto the armrests while getting in and out of the chair.

The tilted chair will catch the weight of the user, allowing them to ease themselves into a natural seated position. This function also allows for a greater sense of independence. The Slimline is rated to 130kg.

Please note: May require assembly.

It is available in three standard sizes:

- Mini to suit height ranges between 122cm (4'10") to 1.57cm (5'2")
- A to suit height ranges between 157cm (5'2") to 173cm (5'8")
- B to suit height ranges between 173cm (5'8") to 180cm (5'11")

The Slimline is proudly made in Australia by Oscar Furniture.

FEATURES

- Modern appearance
- Adjustable headrest
- Adjustable lumbar support
- Customisable
- A large range of modifications available
- Available in 3 base sizes (Mini, A & B)
- Rated to 130kg

WARRANTY

2 Years Warranty: Mechanism, fabric and electronics

5 Years Warranty: Timber Frame and Foam.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au