




OTTOBOCK AVANTGARDE 4 DV WHEELCHAIR

\$5,216.00



SKU: 480F160=0_AA03_V01-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Ottobock Avantgarde 4 DV Wheelchair is a state-of-the-art mobility solution, expertly designed for individuals seeking a high degree of independence and performance. This advanced wheelchair combines cutting-edge technology with user-centric design, offering unparalleled comfort, manoeuvrability, and customisation.

FEATURES

- **Lightweight Frame:** Constructed with a lightweight yet durable material, enhancing ease of movement and transport.
- **Customisable Design:** Offers extensive customisation options in terms of size, colour, and accessories to suit individual preferences and needs.
- **Adjustable Seating System:** Features an adjustable seat for optimal comfort and posture support, accommodating various body types.
- **High Manoeuvrability:** Designed for excellent agility and precision, facilitating smooth navigation in tight spaces.
- **Foldable and Compact:** Easily foldable for convenient storage and transportation, making it ideal for active users.
- **Robust Wheels and Tires:** Comes with high-quality wheels and durable tires suited for various outdoor conditions.
- **Anti-tippers:** The anti-tipper is spring-loaded, with multiple adjustment possibilities. This ensures straightforward operation and the tip-resistance of the wheelchair is assured even with extreme settings.
- **Integrated pockets:** Two pockets below the seat upholstery provide space for the storage of

personal items, protecting them against third-party access.

The Ottobock Avantgarde 4 DV Wheelchair represents a perfect blend of functionality and design, making it a top choice for individuals who value mobility, comfort, and style. Its advanced features and customisation options ensure that it meets the diverse needs and preferences of its users, enhancing their quality of life and independence.

SPECIFICATIONS

Description	Avantgarde 4 DV
SKU	480F160=0_AA03_C
Total Width	490mm - 725mm
Total Length	720mm - 1040mm
Total Height	580mm - 1100mm
Folding Size LxWxH	260mm
Seat Width	40.5cm (15.9"); 46cm (18.1"); 50cm (19.6")
Seat Depth	360mm - 540mm
Seat Height Anterior	380mm - 540mm
Seat Height Posterior	360mm - 510mm
Back Height	250mm - 550mm
Lower Leg Length	160mm - 550mm
Front Angle Frame	75°
Product Weight	9.7kg
SWL	140kg

BROCHURE

[Ottobock Avantgarde 4 Brochure](#)