




## OTTOBOCK AVANTGARDE 4 DV WHEELCHAIR

\$5,216.00 - \$5,668.00



**SKU:** 480F160=0\_AA03\_V01-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The OttoBock Avantgarde 4 DV Wheelchair is a state-of-the-art mobility solution, expertly designed for individuals seeking a high degree of independence and performance. This advanced wheelchair combines cutting-edge technology with user-centric design, offering unparalleled comfort, manoeuvrability, and customisation.

### FEATURES

- **Lightweight Frame:** Constructed with a lightweight yet durable material, enhancing ease of movement and transport.
- **Customisable Design:** Offers extensive customisation options in terms of size, colour, and accessories to suit individual preferences and needs.
- **Adjustable Seating System:** Features an adjustable seat for optimal comfort and posture support, accommodating various body types.
- **High Manoeuvrability:** Designed for excellent agility and precision, facilitating smooth navigation in tight spaces.
- **Foldable and Compact:** Easily foldable for convenient storage and transportation, making it ideal for active users.
- **Robust Wheels and Tires:** Comes with high-quality wheels and durable tires suited for various outdoor conditions.
- **Anti-tippers:** The anti-tipper is spring-loaded, with multiple adjustment possibilities. This ensures straightforward operation and the tip-resistance of the wheelchair is assured even with extreme settings.
- **Integrated pockets:** Two pockets below the seat upholstery provide space for the storage of

personal items, protecting them against third-party access.

The Ottobock Avantgarde 4 DV Wheelchair represents a perfect blend of functionality and design, making it a top choice for individuals who value mobility, comfort, and style. Its advanced features and customisation options ensure that it meets the diverse needs and preferences of its users, enhancing their quality of life and independence.

## SPECIFICATIONS

Description	Avantgarde 4 DV
<b>SKU</b>	480F160=0_AA03_C
<b>Total Width</b>	490mm - 725mm
<b>Total Length</b>	720mm - 1040mm
<b>Total Height</b>	580mm - 1100mm
<b>Folding Size   LxWxH</b>	260mm
<b>Seat Width</b>	40.5cm (15.9"); 46cm (18.1"); 50cm (19.6")
<b>Seat Depth</b>	360mm - 540mm
<b>Seat Height   Anterior</b>	380mm - 540mm
<b>Seat Height   Posterior</b>	360mm - 510mm
<b>Back Height</b>	250mm - 550mm
<b>Lower Leg Length</b>	160mm - 550mm
<b>Front Angle Frame</b>	75°
<b>Product Weight</b>	9.7kg
<b>SWL</b>	140kg

## BROCHURE

[Ottobock Avantgarde 4 Brochure](#)