

## OTTOBOCK LIVO SEAT CUSHION

**\$167.00 - \$205.00** Price range: \$167.00 through \$205.00



**SKU:** 476C00-SK1464-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Ottobock Livo Seat Cushion is a lightweight and durable wheelchair seat cushion that provides excellent pressure redistribution and skin protection for users with moderate to high risk of skin breakdown. It is made with a combination of unique rest suspension foam and ultra-thin LiquiCell packets, which work together to absorb shocks and distribute pressure evenly across the user's seating surface.

The Livo Seat Cushion is also designed to promote airflow and reduce heat buildup, helping to keep the user's skin cool and dry. It has a water-resistant multi-stretch Dartex outer cover that is easy to clean and maintain.

The Livo Seat Cushion is available in a variety of sizes and thicknesses to accommodate different users and needs. It is also compatible with a wide range of wheelchairs and seating systems.

### FEATURES

- Excellent pressure redistribution and skin protection
- Lightweight and durable
- Promotes airflow and reduces heat buildup
- Water-resistant multi-stretch Dartex outer cover
- Available in a variety of sizes and thicknesses
- Compatible with a wide range of wheelchairs and seating systems

The Ottobock Livo Seat Cushion is a great choice for wheelchair users who are looking for a comfortable and supportive seat cushion that can help prevent skin breakdown. It is ideal for users

who spend extended periods of time sitting in their wheelchairs, such as those with spinal cord injuries, multiple sclerosis, or cerebral palsy.

## **SPECIFICATIONS**

- **Materials:** Rest suspension foam and ultra-thin LiquiCell packets
- **Cover:** Water-resistant multi stretch Dartex
- **Weight capacity:** 150kg