




## OTTOBOCK MOTUS XXL WHEELCHAIR

**\$3,929.00 - \$4,606.99** Price range: \$3,929.00 through \$4,606.99



**SKU:** 480F61=2\_AA05\_V01-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Ottobock Motus XXL Wheelchair stands out in the realm of adaptive wheelchairs. It's not just about mobility; it's about delivering a reliable and efficient experience, built around the user's specific needs. With the ability to customise according to preferences and requirements, the Motus XXL promises optimal comfort and enhanced mobility for its users.

### FEATURES

- **Robust Folding Mechanism:** Its unique internal cross-brace folding mechanism is both reliable and durable. This ensures ease of storage and transport, making it perfect for those on the go.
- **Enhanced Weight Support:** With the added option of a double cross brace folding mechanism, the wheelchair can comfortably support a user weight of up to 140kg.
- **Constructed with Excellence:** Made from lightweight aluminium, the chair promises sturdiness without compromising on manoeuvrability.
- **Protection & Comfort:** Integrated side panels protect the user's clothing from any contact with the wheels. Moreover, the breathable backrest upholstery ensures the user remains comfortable, even during prolonged periods of sitting.
- **Ease of Transfer:** Swing-away footplates, available as an option, enable effortless transfers and improved access in front of the chair. It also supports foot propulsion, adding an extra layer of mobility.
- **Travel-friendly Design:** The optional half-folding seat back feature allows the chair to fold forward, reducing its overall height. This makes the Motus XXL an excellent choice for travellers, ensuring it can be stored compactly.

## SPECIFICATIONS

|                              |                                                                                                                |
|------------------------------|----------------------------------------------------------------------------------------------------------------|
| Frame material               | Aluminium                                                                                                      |
| Weight from                  | approx. 13kg (CV)                                                                                              |
| Maximum load capacity        | 125kg (up to 48cm seat width)<br>140kg (from 50.5cm seat width)                                                |
| Centre of gravity adjustment | 48 positions through the rear wheel adaptor                                                                    |
| Wheel locks/brakes           | Knee-level wheel lock<br>Drum brake                                                                            |
| Seat width                   | 45.5cm (17.9"); 50.5cm (19.8")                                                                                 |
| Seat depth                   | 36 - 54cm                                                                                                      |
| Overall length               | 83 - 111cm                                                                                                     |
| Overall width                | 52 - 72cm with standard rear wheels and 0° camber<br>55.5 - 75.5cm with rear wheels with drum brake, 0° camber |
| Back                         | Straight, height adjustable / 30° angle adjustable                                                             |
| Back height                  | 30 - 50cm                                                                                                      |
| Lower leg length             | 34 - 52cm                                                                                                      |