




OTTOBOCK START M2S LIGHTWEIGHT WHEELCHAIR



SKU: 480F74=30102_CT-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Experience mobility like never before with the Ottobock Start M2S Lightweight Wheelchair. As an entry-level model in the renowned Start series, this wheelchair might be your introduction to the Ottobock world, but its superior design and functionality speak volumes about the quality and innovation embedded within it.

FEATURES

- **Sophisticated Design:** Crafted without the use of steel components and sporting a frame identical to other third-generation Start models, its aesthetic is both modern and unobtrusive. The silver metallic frame colour paired with black pads adds a touch of sophistication to its look.
- **High Functionality:** The Ottobock Start M2S doesn't just focus on design but emphasises functionality. The integrated tip assists in making it easier to manoeuvre stairs and provides maximum security even at a basic level.
- **Adjustability and Comfort:** The wheelchair comes equipped with split, pivot-away footrests and side panels featuring height-adjustable armrests. This ensures optimal support and comfort for the user. Moreover, the individual, angle-adjustable lift-up footrests ensure the user can set their preferred foot position.
- **Easy Transport:** With its ability to be folded to a narrow measure in just a few steps, the Ottobock Start M2S is the perfect companion for travel. The wheelchair's hook and loop tape design ensures it remains secure during transport.
- **Quick Customisation:** Tailoring the wheelchair to meet individual requirements is a breeze. While the back height is set at a fixed 42.5 cm, the seat depth, seat height, and other configurations can be quickly adjusted by a rehab technician.

- **Efficient Wheels and Brakes:** The patented wheel lock mechanism requires minimal effort, ensuring safety. Furthermore, rear wheels equipped with quick-release axles facilitate easy and space-efficient loading.
- **Frame and Upholstery:** The wheelchair boasts an aluminium tube frame, which guarantees durability and lightweight manoeuvrability. The seat and back upholstery are crafted from durable black nylon, ensuring long-term use without wear and tear.

The Ottobock Start M2S Lightweight Wheelchair is the epitome of durability, comfort, and efficiency. It effortlessly combines user-friendly features with a sleek design, ensuring you have a reliable mobility companion. Whether you're venturing out into the city or simply manoeuvring around your home, the Ottobock Start M2S is designed to enhance your mobility experience. Step into a world of unparalleled mobility and comfort with the Ottobock Start M2S today.

SPECIFICATIONS

- Seat width:
 - 40.5cm (15.9")
 - 43cm (16.9")
 - 45.5cm (17.9")
 - 48cm (18.8")
 - 50.5cm (19.8")
 - 53.5cm (21")
- Seat depth: Preset to 43cm, can also be set to 40cm or 46cm
- Back height:
 - Standard back: 42.5cm
 - Angle-adjustable back: 47.5cm
- Front seat height: 49cm
- Seat inclination: 3 degrees
- Lower leg length: 31cm - 47cm
- Weight from: 15kg
- Maximum user weight: 125kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA
5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave,
Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA
5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au