

## OTTOBOCK TERRA CUSHION

\$353.00 - \$420.00



SKU: 476C00=SK602-H



## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

When it comes to creating a superior seating experience, the OttoBock Terra Cushion leads the pack. Leveraging internal expertise gained from the automotive seating industry, Ottobock has pioneered a unique foam that ensures unmatched comfort and durability for wheelchair users.

### FEATURES

- **Highly Engineered Foams:** Benefitting from extensive research and experience in automotive seating, Ottobock has crafted a foam that's perfect for prolonged sitting.
- **Innovative Combination:** Merging the flow properties of viscoelastic foam with the resilience of top-tier urethane foams, this cushion introduces a whole new seating experience for wheelchair users.
- **Advanced Molecular Structure:** Distinct from traditional foams, the molecules in Ottobock's foams are elongated with fewer branches. This design promotes enhanced mobility of molecules under various loads.
- **Optimal Immersion & Envelopment:** The cushion guarantees superior pressure redistribution. Its ability to optimize immersion and envelopment ensures utmost comfort.
- **Effective Shock Absorption:** As users navigate different terrains, the unique molecular flow of the foam helps in effectively absorbing impacts and shocks. This provides a spring-like or suspension sensation, greatly enhancing the ride's comfort.
- **Durability:** A standout feature of the OttoBock Terra Cushion is its long-lasting nature. The foam resists breakdown and compression set, ensuring that it retains its shape and comfort even with extended usage.
- **Lightweight Comfort:** Experience the sensation of a suspension without adding any extra

weight to the wheelchair. The cushion provides a spring/suspension feel without the need for additional suspension components.

### **Why Choose the OttoBock Terra Cushion?**

The OttoBock Terra Cushion is not just another wheelchair cushion. It is the outcome of extensive research and innovation. If you are in search of a cushion that offers a seamless blend of comfort, durability, and shock absorption, then look no further. Whether you're navigating the smooth interiors of a building or the rugged paths outdoors, the Terra Cushion promises a comfortable journey. Say goodbye to the discomforts of prolonged sitting and embrace the revolution that is the OttoBock Terra Cushion.

### **SPECIFICATIONS**