



### **DAYS TILTING OVERBED TABLE**

**\$130.00**



**SKU: HLC-751C**

## **PRODUCT DESCRIPTION**

## **PURPOSE AND BENEFITS**

The Days Tilting Overbed Table can be tilted to assist users with seeing and reaching items on the table. This table is height adjustable with a split legged base and has a laminated teak effect top with a rail on the edge to stop items falling off. The swivel castors make moving this overbed table around the bedroom simple and smooth. The frame is chrome plated and powder coated steel.

The table is adjustable in height and comes with swivel castors for ease of manoeuvrability. There are rails on the edges of the tables to prevent items slipping off.

## **SPECIFICATIONS**

- Maximum Load: 10kg
- Height Adjustment: 610mm - 910mm
- Table Size: 610mm x 400mm
- Flatpacked