




## OXFORD ADVANCE PORTABLE FLOOR HOIST



**SKU:** HF2025

## PRODUCT DESCRIPTION

## PURPOSE AND BENEFITS

The Oxford Advance folding mobile patient lift is one of the most functional and stylish in its class today. Especially suited to home care and small nursing institutions where space sometimes is simply not available, the Advance also allows a user to travel outside of the confines of their own home.

## FEATURES

- Oversized handle for increased control
- Patent pending push foot pad reduces the force needed to initiate movement
- Compact design
- Swan neck leg design accommodates large furniture
- Patent pending triangular folded design allows the lift to securely stand unaided
- No tools required to fold
- All aluminium construction
- Foot operated leg spreading mechanism requires no bending to operate

## SPECIFICATION

<b>Length</b>	1300 mm
<b>Maximum Height</b>	1860 mm
<b>Minimum Height</b>	1360 mm
<b>Height (Folded)</b>	450 mm

<b>Depth (Folded)</b>	1180 mm
<b>Width</b>	550 mm
<b>Height at Maximum Reach</b>	1170 mm
<b>Reach at Maximum Height</b>	650 mm
<b>Reach at Minimum Height</b>	350 mm
<b>Turning Radius</b>	1420 mm
<b>Legs Open (External Width)</b>	1080 mm
<b>Legs Open (Internal Width)</b>	1000 mm
<b>Legs Closed (External Width)</b>	670 mm
<b>Legs Closed (Internal Width)</b>	570 mm
<b>Overall Height of Legs</b>	115 mm
<b>Front Castor</b>	75 mm
<b>Rear Braked Castor</b>	100 mm
<b>Maximum Weight Capacity</b>	155 kg
<b>Brand</b>	Oxford

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
 Brighton SA 5048  
 08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
 Ave, Elizabeth Vale SA 5112  
 08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
 Reynella SA 5161  
 08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
 03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
 03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)