

## OXFORD ASCEND STAND-AID LIFT

\$7,174.00



**SKU:** HF3320



## PRODUCT DESCRIPTION

### Purpose & Benefits of Oxford Ascend Stand-Aid Lift

The **Oxford Ascend Stand-Aid Lift** is a compact and versatile patient lifter designed to support safe and efficient standing and transfer movements for individuals with limited mobility. Ideal for both home and facility use, the Ascend combines lift and transfer functionality in a single, easy-to-manooeuvre unit. Its ergonomic design allows caregivers to safely assist users in standing, pivoting, or moving from bed to chair, commode, or other seated surfaces with minimal strain.

With an adjustable knee pad, padded leg support, and powered lift mechanism, the Ascend promotes user comfort while reducing the risk of injury for both clients and caregivers. It is well-suited for individuals who have some weight-bearing ability and can actively participate in transfers, making it a practical and dignified alternative to full hoists.

### Key Features:

- Powered stand-assist lift for active transfers
- Compact design for easy use in tight spaces

- Adjustable knee pad and padded leg support for added comfort
- Ergonomic handle and push bar for improved caregiver control
- Compatible with Oxford standing slings
- Safe working load of up to 170 kg
- Easy-to-use handset and emergency lowering function
- Ideal for home care, hospitals, and aged care facilities

**Specifications:**

- Safe working load: 200 kg
- Maximum overall length: 1100 mm
- Minimum overall length: 1000 mm
- Maximum overall height: 1710 mm
- Minimum overall height: 1070 mm
- Maximum height to attachment point: 1665 mm
- Minimum height to attachment point: 990 mm
- Turning radius: 1320 mm
- Legs open - external width: 995 mm
- Legs open - internal width: 885 mm
- Legs closed - external width: 630 mm
- Legs closed - internal width: 510 mm
- Overall height of legs: 120 mm
- Ground clearance: 35 mm
- Front twin castors: 100 mm
- Rear braked castors: 100 mm