



## **OXFORD JOURNEY FOLDING STAND HOIST**

**\$6,179.00**



**SKU: HF2027**

## **PRODUCT DESCRIPTION**

### **PURPOSE AND BENEFITS**

The Oxford Journey lift is a compact foldable stand with high functionality. Its sleek design allows for the Oxford lift to be compact and lightweight which makes it easy for carers to move the lift through tight spaces and allows for easy storage.

The first on the market adjustable sling attachment cradle mechanism supports a range of patient heights and sizes, allowing your facility to have one lift that meets the needs of many. This careful attention to product design with the Journey ensures that both the caregiver's and the resident's needs have been met.

### **FEATURES**

- Compact folding design for easy storage and transportation
- Supports a range of patient heights and sizes
- Unique adjustable cow-horn function
- Electric hand control operation
- Emergency operation
- Smart™ Monitor diagnostic control system
- Contoured, height adjustable knee pad
- Removable foot tray for walking|rehabilitation activity
- Easily manoeuvrable with 125mm rear locking castors and 75mm front twin swivel castors
- Sling NOT included with purchase

### **SPECIFICATIONS**

|                              |               |
|------------------------------|---------------|
| <b>Length</b>                | 1000mm        |
| <b>Width</b>                 | 685mm         |
| <b>Height</b>                | 1130mm-1600mm |
| <b>Leg Height</b>            | 109mm         |
| <b>Legs Open   Inside</b>    | 876mm         |
| <b>Legs Closed   Outside</b> | 689mm         |
| <b>Lift Range</b>            | 7200mm-1530mm |
| <b>Leg Gait Range</b>        | 590mm-875mm   |
| <b>Electric Leg</b>          | True          |
| <b>Weight</b>                | 40.2kg        |
| <b>SWL</b>                   | 155kg         |