



OXFORD JOURNEY FOLDING STAND HOIST



SKU: HF2027

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Oxford Journey lift is a compact foldable stand with high functionality. Its sleek design allows for the Oxford lift to be compact and lightweight which makes it easy for carers to move the lift through tight spaces and allows for easy storage.

The first on the market adjustable sling attachment cradle mechanism supports a range of patient heights and sizes, allowing your facility to have one lift that meets the needs of many. This careful attention to product design with the Journey ensures that both the caregiver's and the resident's needs have been met.

FEATURES

- Compact folding design for easy storage and transportation
- Supports a range of patient heights and sizes
- Unique adjustable cow-horn function
- Electric hand control operation
- Emergency operation
- Smart™ Monitor diagnostic control system
- Contoured, height adjustable knee pad
- Removable foot tray for walking|rehabilitation activity
- Easily manoeuvrable with 125mm rear locking castors and 75mm front twin swivel castors
- Sling NOT included with purchase

SPECIFICATIONS

Length	1000mm
Width	685mm
Height	1130mm-1600mm
Leg Height	109mm
Legs Open Inside	876mm
Legs Closed Outside	689mm
Lift Range	7200mm-1530mm
Leg Gait Range	590mm-875mm
Electric Leg	True
Weight	40.2kg
SWL	155kg

For further information go to [medimart.com.au](https://www.medimart.com.au) or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
 1300 633 468

support@medimarnt.com.au

