




OXFORD MIDI ELECTRIC LIFTER



SKU: HF2029

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Oxford Midi 180 has become the 'gold standard' in homecare, offering great performance for a hoist of its size. It's the perfect solution for both community and nursing environments while still being sympathetic to space constraints it might encounter.

FEATURES

- 180kg Safe Working Load
- Smart Monitor diagnostic and service controls as standard
- Foot pedal leg opening mechanism with latch/lock functionality
- Electric leg option available
- High-grade Tente castors for easy manoeuvrability
- Detachable battery pack and desktop charger
- Optional Class III digital weigh-scale

SPECIFICATIONS

Safe Working Load: 180kg

Length: 1225mm

Min Overall Height: 1240mm

Max Overall Height: 1840mm

Min Spreader Bar Height: 525mm

Max Spreader Bar Height: 1660mm

Spreader Bar Reach at Min Height: 440mm

Spreader Bar Reach at Max Height: 420mm

Spreader Bar Height at Max Reach: 1060mm

Max Reach: 656mm

Turning Radius: 1235mm

Legs Open - External Width: 1170mm

Legs Open - Internal Width: 1070mm

Legs Closed - External Width: 600mm

Legs Closed - Internal Width: 540mm

Overall Height of Legs: 100mm

Ground Clearance: 25mm

Front Twin Castors: 75mm

Rear Braked Castors: 100mm

Total Weight: 42kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au