




OXFORD MIDI ELECTRIC LIFTER



SKU: HF2029

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Oxford Midi 180 has become the 'gold standard' in homecare, offering great performance for a hoist of its size. It's the perfect solution for both community and nursing environments while still being sympathetic to space constraints it might encounter.

FEATURES

- 180kg Safe Working Load
- Smart Monitor diagnostic and service controls as standard
- Foot pedal leg opening mechanism with latch/lock functionality
- Electric leg option available
- High-grade Tente castors for easy manoeuvrability
- Detachable battery pack and desktop charger
- Optional Class III digital weigh-scale

SPECIFICATIONS

Safe Working Load: 180kg

Length: 1225mm

Min Overall Height: 1240mm

Max Overall Height: 1840mm

Min Spreader Bar Height: 525mm

Max Spreader Bar Height: 1660mm

Spreader Bar Reach at Min Height: 440mm

Spreader Bar Reach at Max Height: 420mm

Spreader Bar Height at Max Reach: 1060mm

Max Reach: 656mm

Turning Radius: 1235mm

Legs Open - External Width: 1170mm

Legs Open - Internal Width: 1070mm

Legs Closed - External Width: 600mm

Legs Closed - Internal Width: 540mm

Overall Height of Legs: 100mm

Ground Clearance: 25mm

Front Twin Castors: 75mm

Rear Braked Castors: 100mm

Total Weight: 42kg