



OXFORD PRESENCE FLOOR LIFTER



SKU: LM-3499

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Oxford Presence is a heavy-duty mobile patient lift, competent enough to handle almost any patient-handling task. Its incredible lifting range ensures it can cope easily with the most critical of transfers, from the floor to high level surfaces.

FEATURES

- Safe working load 227kg
- Outstanding lifting range (440mm – 1910mm)
- SmartTM Monitor diagnostic and service controls as standard
- Interchangeable spreader bar and cradle options
- Ergonomic design considerations and interaction points
- Powered leg positioning for excellent access around furniture
- Low leg model available especially suited to low bed frames
- Powered cradle and digital weigh-scale options 5 yr warranty on aluminium Frame. 2 yr warranty on mechanical and electrical components. Battery and casters 90 day warranty

SPECIFICATIONS

SAFE WORKING LOAD	227 kgs
MAXIMUM OVERALL LENGTH	1460 mm
MINIMUM OVERALL LENGTH	1400 mm
MAXIMUM OVERALL HEIGHT	2065 mm
MINIMUM OVERALL HEIGHT	1460 mm

SPREADER BAR MAXIMUM HEIGHT	1910 mm
SPREADER BAR MINIMUM HEIGHT	430 mm
HEIGHT AT MAXIMUM REACH	1310 mm
REACH AT MAXIMUM HEIGHT	650 mm
REACH AT MINIMUM HEIGHT	400 mm
MAXIMUM REACH	930 mm
TURNING RADIUS	1630 mm
LEGS OPEN - INTERNAL WIDTH	1020 mm
LEGS CLOSED - EXTERNAL WIDTH	670 mm
FRONT TWIN CASTERS	100 mm
REAR BRAKED CASTERS	100 mm
OVERALL HEIGHT OF LEGS	120 mm
WITH LOW LEG OPTION	60 mm
GROUND CLEARANCE	35 mm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au