



OXFORD PRESENCE FLOOR LIFTER

\$6,332.00



SKU: HF2031

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Oxford Presence is a heavy-duty mobile patient lift, competent enough to handle almost any patient-handling task. Its incredible lifting range ensures it can cope easily with the most critical of transfers, from the floor to high level surfaces.

FEATURES

- Safe working load 227kg
- Outstanding lifting range (440mm – 1910mm)
- SmartTM Monitor diagnostic and service controls as standard
- Interchangeable spreader bar and cradle options
- Ergonomic design considerations and interaction points
- Powered leg positioning for excellent access around furniture
- Low leg model available especially suited to low bed frames
- Powered cradle and digital weigh-scale options 5 yr warranty on aluminium Frame. 2 yr warranty on mechanical and electrical components. Battery and casters 90 day warranty

SPECIFICATIONS

SAFE WORKING LOAD	227 kgs
MAXIMUM OVERALL LENGTH	1460 mm
MINIMUM OVERALL LENGTH	1400 mm
MAXIMUM OVERALL HEIGHT	2065 mm
MINIMUM OVERALL HEIGHT	1460 mm

SPREADER BAR MAXIMUM HEIGHT	1910 mm
SPREADER BAR MINIMUM HEIGHT	430 mm
HEIGHT AT MAXIMUM REACH	1310 mm
REACH AT MAXIMUM HEIGHT	650 mm
REACH AT MINIMUM HEIGHT	400 mm
MAXIMUM REACH	930 mm
TURNING RADIUS	1630 mm
LEGS OPEN - INTERNAL WIDTH	1020 mm
LEGS CLOSED - EXTERNAL WIDTH	670 mm
FRONT TWIN CASTERS	100 mm
REAR BRAKED CASTERS	100 mm
OVERALL HEIGHT OF LEGS	120 mm
WITH LOW LEG OPTION	60 mm
GROUND CLEARANCE	35 mm