


## P5500 PULL-APART BED



SKU: P5500-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Caremed Alrick P5500 Pull-Apart Bed is a highly functional and versatile medical bed designed to cater to the needs of patients and healthcare providers in various medical settings. This bed is specially crafted to offer ease of use, comfort, and practicality, particularly in environments where space optimization and patient mobility are key considerations.

### FEATURES

- **Pull-Apart Functionality:** The unique pull-apart feature allows the bed to be easily separated into two parts, facilitating efficient transport and storage. This is particularly useful in settings where space is limited or the bed needs to be moved frequently.
- **Adjustable Height:** Equipped with an electric height adjustment mechanism, the bed can be raised or lowered to suit the specific needs of patients and caregivers, enhancing comfort and facilitating easier patient transfers.
- **Robust Construction:** Built with durable materials, the P5500 is designed to withstand regular use in a healthcare environment, ensuring longevity and reliability.
- **User-Friendly Controls:** The bed comes with intuitive controls that make adjustments simple and efficient, reducing the workload for caregivers and improving the overall patient experience.
- **Safety Features:** Includes essential safety elements like side rails to prevent patient falls and lockable casters to secure the bed in place during use.
- **Comfortable Sleeping Surface:** The bed is designed with a comfortable mattress support platform, suitable for various types of medical mattresses, ensuring patient comfort during recovery or treatment.

The Caremed Alrick P5500 Pull-Apart Bed is an excellent choice for medical facilities looking for a bed that offers flexibility, comfort, and ease of use. Its innovative design and robust construction make it a valuable asset in providing high-quality care and enhancing the efficiency of healthcare operations.

## **SPECIFICATIONS**

- **Size Specifications:** King Single
- Overall Width - 1070mm
- Overall Length - 2190mm
- Sleeping Surface - 1050 x 2000mm
- Height - 295-740mm
- Safe Working Load - 200kg
  
- **Size Specifications:** Single
- Overall Width - 930mm
- Overall Length - 2190mm
- Sleeping Surface - 900 x 2000mm
- Height - 295-740mm
- Safe Working Load - 200kg